


ML 611-681

ML 681

Dati tecnici

Technical data

ML611 - ML681

Diametro flacone mm	35-90	Bottle diameter mm
Altezza flacone mm	65-200	Bottle height mm
Velocità meccanica max (cicli/min.)	20-100	Max mechanical speed (cycles/minute.)

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Linear Filler ML 611 & Capper-Pluggger ML 681

Linear Filling Machine ML 611

Automatic intermittent motion machine suitable to fill bottles with volumetric dosing through ELECTROMAGNETIC FLOW METERS or CORIOLIS MASS FLOW METERS. Cantilevered (balcony type) machine layout to better comply with CGMP directives. Clear separation between the bottle transport and mechanical area. Bottle transport on structure made of stainless steel. Machine movements achieved by electric motor and Inverter with speed settings from Display. One first infeed tabletop conveyor and a second discharge tabletop conveyor. Rapid size changes over with a minimum Number of parts to replace or to position.

Filling Systems: Various filling circuits with different types of measuring equipment are available.

- ELECTROMAGNETIC FLOW METERS, (flow meter based on electromagnetic conductivity of products) or
- CORIOLIS MASS FLOW METERS (mass flow meters based on the Coriolis principle).

The number of dosing units, from 1 to 10 is dedicated to the volume and speed requirements. The bottle transfer from the first infeed conveyor onto the elevating platform under the filling station is driven by a brushless motor with variable stroke, depending on the size to be handled. Size configuration is memorized on the Computer.

The variable nozzle stroke is commanded by a brushless motor and can be set on the PC for best bottom-up filling sequence. Filling circuits are designed to accept a CIP-SIP (Clean in Place or Steam in Place) system as Optional. Product contact parts as well as the nozzle holding rake are of easy replacement.

Operator PC Color Display Panel for:

- Machine function management.
- Operator interface.
- Size configuration and memorization.
- Dosing volume settings.
- Automatic adjustment of nozzle height.
- Programmed maintenance.

Linear, Robotized Capping and Plugging Machine ML 681

The ML 681 has been designed to comply with the possibility to perform the assembling of different closing components as well as their inserting, capping and plugging function. Typical applications relate to spray pumps, screw caps and cover (dust) caps. The characteristics point in particular to the Cosmetic Industry, which utilize components of dedicated form and application and may not be handled efficiently with traditional closing machines. Cantilever (balcony) design corresponding to CGMP directives. The machine may be equipped with up to THREE robotized units, depending on the typology of the closing components. The container transport is achieved by tabletop conveyors with independent motor and Inverter. A screw synchronizes container infeed and discharge.

Operator PC Color Display Panel for:

- Machine function management.
- Operator interface.
- Size configuration and memorization.
- Programmed maintenance.

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Dosing group with Coriolis Mass Flow Meters. Product mass flow meters suitable for liquids without electromagnetic conductivity. Nozzle holding rake with bottle neck centering cones.

Picking of spray pump, straightening-centering of the pump under the screw capping head. Both movements are commanded by brushless motor, all size parameters (stroke, capping speed etc.) are managed by the PC.

Robotized final screw capping with torque control from Display as well as programmable stroke adequate to the pump length.

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5

6

- Robotized head with two movements. The first picks the press-on cap and places it onto the cover cap. The second places the cover cap onto the container. (See picture 5). Adjustable stroke based on the height of the caps. Cap feeding conveyors driven by brushless motors.*

- The assembled cover, with a press-on movement, is placed onto the container. A control verifies the correct positioning, otherwise triggers the reject.*

- View of the cover feeding unit, where the cover is taken off the chute. The vibratory is placed at a low height of 1150mm from floor level.*